

## GeneProof Factor V Leiden PCR Kit



*In vitro* diagnostic medical device

The kit has been manufactured according to EC Directive 98/79/EC as an *in vitro* diagnostic medical device and it has been designed for professional use in specialized clinical and research laboratories.

### KIT CONTENT

REF	FV/050 50 rxn	FV/100 100 rxn
<b>MasterMix</b> FV Leiden	2x450 µl	4x450 µl
<b>Positive Control</b> FV Leiden wild type	1x50 µl	1x50 µl
<b>Positive Control</b> FV Leiden mutant	1x50 µl	1x50 µl
<b>Positive Control</b> FV Leiden heterozygote	1x50 µl	1x50 µl

### STORAGE AND TRANSPORTATION CONDITIONS

The kits should be transported and stored at temperatures between -85 °C and -10 °C. The kit will remain stable at least until the expiry date printed on the package, if the storage temperature is kept. Repeated freezing and thawing of the kit components may result in lower detection quality.

### TECHNICAL SPECIFICATION

<b>Target sequence</b>	DNA sequence of a gene encoding hemocoagulation factor V (Leiden variation)
<b>Specificity</b>	G1691A single-nucleotide mutation in exon 10 (Arg 506 Gln)
<b>Qualitative detection</b>	mutation present x not present
<b>Evaluation</b>	G/G - standard homozygote (wild type) A/A - mutant homozygote G/A - heterozygote
<b>Sample types</b>	EDTA whole blood
<b>Quality Control</b>	regularly tested by Instand e.V. External Quality Assessment Panel

## METHOD PRINCIPLES

The PCR kit is designed to detect G1691A mutation in gene for human factor V Leiden by the real-time polymerase Chain Reaction (PCR) method. Method is based on amplification and detection of target sequence using allele specific fluorophore labelled probes. Target sequence is single nucleotide replacement of adenine for guanine in site 1691 (G1691A). Presence of wild-type allele (G1691G) is indicated in FAM fluorescent channel and mutant allele (A1691A) in HEX fluorescent channel. In case of heterozygous genotype (G1691A) signal is detected in both channels. Detection kit contains Ready to Use MasterMix and utilizes "hot start" technology, minimizing non-specific reactions and assuring maximum sensitivity. Kit is designed for *in vitro* diagnostic.

## GENETIC DIAGNOSTIC TECHNOLOGY


### / SAMPLE

sample contains  
— wild-type allele  
— mutant allele





### / ISOLATION

the extracted sample DNA is added into the Ready to Use MasterMix and the tube is inserted into the real-time device



### / PCR AMPLIFICATION

individual alleles are amplified in the course of PCR



### / EVALUATION

#### WILD-TYPE HOMOZYGOTE

– only a wild-type allele in the FAM channel is detected

#### HETEROZYGOTE

– both alleles are detected:

1. a wild-type allele in the FAM channel
2. a mutant allele in the HEX channel

#### MUTANT HOMOZYGOTE

– only a mutant allele in the HEX channel is detected



## USER MANUAL

### SAMPLING AND SAMPLE STORAGE

Sampling does not require any special preparation of the patient. Usually 3 ml of peripheral blood are sampled into a tube with EDTA. It is possible to keep the blood in a refrigerator at the temperature between +2 °C and +8 °C for 3 days; for long-term storage keep the samples at the temperature between -85 °C and -10 °C.

### NUCLEIC ACID PURIFICATION

Nucleic acid isolation should be performed by isolation kits available at the market according to protocols for the particular clinical material isolation. The manufacturer recommends the following isolation kits:

GeneProof PathogenFree DNA Isolation Kit  
croBEE NA16 Nucleid Acid Extraction System

### PCR SETUP

1. Add 18 µl of MasterMix into PCR tubes.
2. Add 2 µl of the DNA sample or 2 µl of Positive Control into the individual PCR tubes. The final reaction mix volume will be 20 µl. *It is necessary to keep all components at +2°C to +8°C during the PCR preparation.*
3. Close the tubes, centrifuge shortly, insert them into the device and let them amplify according to the following PCR profile. *Be very careful when handling the Positive Control or the clinical material, incorrect handling could result in contamination and the consequent impairment of the kit components or the MasterMix! The manufacturer is not responsible for the kit impairment due to incorrect handling.*

### AMPLIFICATION PROFILE

Step	Temperature	Time	Data collection	Cycles
1. Hold	95 °C	10 min		1
	95 °C	10 s		
2. PCR	64 °C	20 s	FAM+HEX	40
	72 °C	20 s		

### VALIDATED INSTRUMENTS









GeneProof PCR kits are designed for use with real-time devices from various manufacturers. This PCR kit has been validated with the following devices:

Applied Biosystems 7300/7500 Real-Time PCR System  
Dx/CFX96™/CFX Connect™ Real-Time PCR Detection System  
LightCycler® 2.0, LightCycler® 480  
LineGene 9600  
Rotor-Gene 3000, Rotor-Gene 6000  
SLAN® Real-Time PCR System  
StepOne™ Real-Time PCR System

GeneProof diagnostic kits are continually validated with various types of devices. Please request the current list at [support@geneproof.com](mailto:support@geneproof.com).



## CLINICAL SAMPLES ANALYSIS EVALUATION

Channel FAM	Channel HEX	Result	Interpretation
		Valid	G/G – standard homozygote (wild type)
		Valid	A/A – mutant homozygote
		Valid	G/A - heterozygote
		Invalid	

### WARNING

A single valid package insert for a specific kit is included in the package or to be requested for the particular lot from the manufacturer. The kit should be disposed of after use according to the current legal regulations considering the fact that the kit doesn't contain any dangerous, infectious or toxic components that would be subject to special safety regulations, and the packaging materials are made of paper and polypropylene. If you have any questions please contact our Customer Service.

#### Customer care and technical support

Tel.: +420543211679  
Fax: +420516770824  
email: support@geneproof.com

#### Orders

Tel.: +420543211679  
Fax: +420516770824  
email: sales@geneproof.com

